

Genotoul Metatoul

ANALYSIS CATALOGUE

Axiom platform



 MetaToul-MetaboHUB
platform



MetaToul is co-funder of
National Infrastructure
MetaboHUB

UNTARGETED METABOLOMICS

Extraction SOP & Automation



Study of Human and Animal Matrices



Analyses
NMR and/or LC-HRMS



Metabolites
Polar and Non-Polar

Internal Databases



Identification
of ~1000 metabolites

Dipeptides	Amino acids
Vitamins	Nucleotides
Catecholamines	Purines
Glutathiones	Pyrimidines
Carnitines	Organic Acids
Steroids	Polyamines
Phospholipids	Glucides
Bile acids	Alcohols
Fatty Acids	
Short-chain metabolites	
Ketones	

Dosage* :
Aldehydes

MS
NMR
MS/NMR

Quantification
in ~60 Metabolic
Subnetworks



*Absolute quantification




Multivariate Statistical
Analyses

Submit a project here




XENOBIOTICS SCREENING



LC-MS or GC-MS analysis of xenobiotics and/or their metabolites following their extraction from biological matrices.

**Relative quantification by HRMS and identification by MS/MS.
Absolute quantification by LC-tQMS in the particular case of
HAAs in meat.**



Mass-Spectrometry Imaging



STATISTICS

Multivariate analysis

- PCA
- PLS-DA
- Repeated measures (A-SCA, AMOPLS)




Single block analysis


Individual analysis of each metabolomics data

Univariate analysis

- Hypothesis testing
- Variable selection (ongoing)



Discrimination of individuals with a metabolic signature



Block contribution to factor prediction

Multi-block analysis

Combination of different analytical metabolomics data to adjust predictive models thanks to a deeper metabolic coverage




DIABLO (mixOmics)
Consensus OPLSDA (Boccard et al.)
Regularized Consensus OPLSDA (ongoing)




The axiom core facility develop and apply analytical chemistry methods by NMR and MS, for analysis of metabolome, xenobiotics and their metabolites, and exposome in toxicology, human and animal health.

Our resources include 17 scientists, 1 NMR, 6 mass spectrometers and 1 robot for sample preparation

2 x LC-MS (tQ)




600 MHz NMR




LC-HRMS (QToF)




GC-MS (Q)



MALDI-LC-IM-HRMS (TIMS-QToF)



LC-HRMS (IQ-X)



Sample preparation robot